FR
DE
EN

Lucidis 108MT - Lucidis 124MT

Instructions d'utilisation
Gebrauchsanweisung
Instructions for Use

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Lucidis Toric intraocular lens

EDOF posterior chamber intraocular lens

(UV-absorbing, foldable hydrophilic acrylic IOL)

Description

Lucidis Toric is a one-piece, flexible intraocular lens implantable in patients in place of the existing crystalline lens which has been clouded by cataracts and affected by corneal astigmatism. The lens is inserted into the capsular bag inside the eye by a surgical procedure after the extraction of the cataract. It consists of a central biconvex optical zone (opt) of 6.0 mm diameter, whose anterior surface has an aspherical structure and two outer parts (haptics) for fixing in the capsular bag. The lens is available with external diameters of 10.8 mm and 12.4 mm. Lucidis Toric lenses are available with a spherical equivalent power of +3.0 to +30.0 diopters and a cylinder power of 1.0 to 4.5 diopters. This information is printed on the product label. The optical equivalent power corresponds to the average power of the lens for distant vision and the cylinder power corresponds to the difference of power between the flat and steep meridian of the lens. The additional depth of focus is designated by «EDOF» and expresses the additional power for near and intermediate vision. The lens is made of a hydrophilic polymer material. The polymer absorbs ultraviolet radiation with a transmittance of less than 10 % for wavelengths under 370 nm (see Fig. 1).

Directions

Lucidis Toric lenses are prescribed for patients suffering from cataracts and pre-existing corneal astigmatism and are requiring surgery. The required powers must be within the range of the diopters available in the Lucidis Toric product portfolio.

Contraindications

Prior to the implantation of an Lucidis Toric lens, a pre-operative assessment should be carried out by a surgeon in order to consider the risk/benefit ratio for the following conditions:

- When the patient already has a disease or physiological conditions that may be exacerbated by the introduction of the lens, or where the lens may interfere with the possibility of examination or treatment of diseases.
- Acute eye disease or external or internal infection
- Fuchs' proliferative diabetic retinopathy
- Severe corneal dystrophy
- Severe myopia
- Chronic or recurrent severe uveitis
- Severe optic nerve atrophy
- Congenital bilateral cataract or non-age related cataract
- Glaucoma
- Choroidal hemorrhage
- Chronic or severe vitreous detachment
- Irregular corneal astigmatism
- Extremely shallow anterior chamber
- Severe micro- or macrophthalmos
- Aniridia
- Retinal detachment
- Macular degeneration
- Microcupill

- Remove the lens from the packaging with a pair of forceps. Do not reverse the lens; the surface visible on opening the blister pack is the anterior surface.
- Immediately position the lens in the middle of the cartridge, with the haptic near the channel of the cartridge pointing left (anterior surface of the lens upward).
- Use open forceps in order to completely avoid touching the optical surface of the lens. Gently fold both wings of the cartridge by pressing on the edges of the optical surface. When the cartridge is nearly closed, pull the forceps back, maintaining only slight pressure to prevent the lens from moving. Make sure that the haptics are properly inserted into the cartridge. If necessary, gently press them using the jaws of the forceps.
- Keep the two wings of the cartridge closed.
- Insert the cartridge into the injector body. Carefully move the plunger of the injector and gently push the lens all the way along the channel. The orientation of the injector tip towards the eye, the pressure applied to the plunger and the speed of injection are all at the discretion of the surgeon.

For instructions on loading the lens into the injector and implementation, refer to the instruction sheet provided with the injector. Recommended injection speed and corresponding incision size (the tests were conducted with the VISCOJECT® 2.2 eye system and METHYLVISUC® solution). Check that the orientation of the Lucidis Toric lens is correct and that it has not been reversed. The anterior surface of the lens should be directed toward the cornea (see Fig. 2). From this point, the operation may continue according to standard surgical protocols. Post-operative follow-up of the patient can also follow stand- ard surgical protocols.

Precautions

- Prior to implantation, check that the model and the power indicated on the label of the lens packaging are correct.
- Do not use the lens if the packaging is damaged or if the expiry date has passed. The sterility of the product may be compromised.
- Handle the lens by the haptics. Avoid all contact with the optical parts.
- SAV-IOL recommends only using injectors designed specifically for the insertion of intraocular lenses.

- The position and size of the incision influences the value of astigmatism and must be considered by the surgeon.
- A misalignment of the intended axis of placement might occur if there is residual viscoelastic or an over-inflation of the capsular bag.
- Severe complications during surgery
 - Choroidal hemorrhage
 - Vitreous loss

Complications

Adverse reactions and complications due to or following surgery may include, but are not limited to:

- Endophthalmitis
- Discharge from the wound, inadequate healing of incision, iris prolapse
- Anterior chamber hemorrhage
- Choroidal hemorrhage
- Hyphema
- Vitreous loss
- Traumatic iritis
- Rupture of the posterior capsule/zone
- Rupture of the haptics during insertion or in the eye
- Difficulties folding the implant/damage to the implant during implantation

Handling precautions

- Do not re-sterilize the lens.
- The lenses should not be exposed to any solutions other than those normally used for standard irrigation (e.g. isotonic saline solution, BSS, viscoelastic, etc.).
- The lens should be handled carefully at room temperature. No attempts should be made to reshape or cut any portion of the lens or to apply undue pressure to the optical portion of the lens with a sharp object.
- Do not allow the lens to dry in air. The lens should be stored in a sterile BSS solution during surgery.

- Lucidis Toric lens must be implanted by an ophthalmic surgeon.
- Store at room temperature, 20° not store in direct sunlight, do not expose to temperatures over 45 °C or under 5 °C for prolonged periods.
- Contamination of the lens may lead to injury or illness of the patient.

Calculation of the power of the Lucidis Toric lens

Prior to the operation to implant the lens, the dioptric spherical equivalent power to be implanted must be determined. Emmetropia is not necessarily a desirable postoperative goal, and factors such as the visual status of the other eye and the life-style of the patient must be considered when determining the lens power to be used. The nominal A-constant is indicated on the product label. Preoperative calculation of the required spherical equivalent lens power should be determined by the surgeon, based on clinical experience and intraoperative placement. The nominal A-constant determines the spherical equivalent and cylindrical power of the lens to be implanted.

Placement of Lucidis Toric lens

The location and size of the incision may affect the axis and/or postoperative corneal astigmatism. To select the appropriate IOL, optical cylinder power, and axis placement, SAV-IOL provides an online Toric Calculator available at www.sav-iol.com. For optimal results, surgeons should customize their surgically induced corneal astigmatism values based upon past results and individual surgical techniques.

The surface of the lens is marked with two indentations that identify the flat meridian to ensure the correct placement of the lens within the capsular bag. The axis marks must be aligned with the intended axis of placement.

Lens packaging

Each Lucidis Toric lens is provided in a double sterile barrier packaging, in a primary packaging containing the BSS, integrated in a pocket and placed in a box with labels and the product information. This device was steam sterilized.

Administration and instructions for use

Lenses must only be implanted by a surgeon who has the skills required for this surgical technique. The standard implantation procedure can be used. General instructions for the use of an injector are included in this document for information purposes.

- Remove the lens from the packaging with a pair of forceps. Do not reverse the lens; the surface visible on opening the blister pack is the anterior surface.
- Immediately position the lens in the middle of the cartridge, with the haptic near the channel of the cartridge pointing left (anterior surface of the lens upward).
- Use open forceps in order to completely avoid touching the optical surface of the lens. Gently fold both wings of the cartridge by pressing on the edges of the optical surface. When the cartridge is nearly closed, pull the forceps back, maintaining only slight pressure to prevent the lens from moving. Make sure that the haptics are properly inserted into the cartridge. If necessary, gently press them using the jaws of the forceps.
- Keep the two wings of the cartridge closed.
- Insert the cartridge into the injector body. Carefully move the plunger of the injector and gently push the lens all the way along the channel. The orientation of the injector tip towards the eye, the pressure applied to the plunger and the speed of injection are all at the discretion of the surgeon.

For instructions on loading the lens into the injector and implementation, refer to the instruction sheet provided with the injector. Recommended injection speed and corresponding incision size (the tests were conducted with the VISCOJECT® 2.2 eye system and METHYLVISUC® solution). Check that the orientation of the Lucidis Toric lens is correct and that it has not been reversed. The anterior surface of the lens should be directed toward the cornea (see Fig. 2). From this point, the operation may continue according to standard surgical protocols. Post-operative follow-up of the patient can also follow stand- ard surgical protocols.

Precautions

- Prior to implantation, check that the model and the power indicated on the label of the lens packaging are correct.
- Do not use the lens if the packaging is damaged or if the expiry date has passed. The sterility of the product may be compromised.
- Handle the lens by the haptics. Avoid all contact with the optical parts.
- SAV-IOL recommends only using injectors designed specifically for the insertion of intraocular lenses.

- The position and size of the incision influences the value of astigmatism and must be considered by the surgeon.
- A misalignment of the intended axis of placement might occur if there is residual viscoelastic or an over-inflation of the capsular bag.
- Severe complications during surgery
 - Choroidal hemorrhage
 - Vitreous loss

Complications

Adverse reactions and complications due to or following surgery may include, but are not limited to:

- Endophthalmitis
- Discharge from the wound, inadequate healing of incision, iris prolapse
- Anterior chamber hemorrhage
- Choroidal hemorrhage
- Hyphema
- Vitreous loss
- Traumatic iritis
- Rupture of the posterior capsule/zone
- Rupture of the haptics during insertion or in the eye
- Difficulties folding the implant/damage to the implant during implantation

Handling precautions

- Do not re-sterilize the lens.
- The lenses should not be exposed to any solutions other than those normally used for standard irrigation (e.g. isotonic saline solution, BSS, viscoelastic, etc.).
- The lens should be handled carefully at room temperature. No attempts should be made to reshape or cut any portion of the lens or to apply undue pressure to the optical portion of the lens with a sharp object.
- Do not allow the lens to dry in air. The lens should be stored in a sterile BSS solution during surgery.

- Lucidis Toric lens must be implanted by an ophthalmic surgeon.
- Store at room temperature, 20° not store in direct sunlight, do not expose to temperatures over 45 °C or under 5 °C for prolonged periods.
- Contamination of the lens may lead to injury or illness of the patient.

SAV-IOL SA
Route des Falaises 74
2000 Neuchâtel
Svizzera
info@sav-iol.com

Main telephone number: +41 32 566 54 00
Regulatory affairs telephone number: +41 32 566 54 08
Regulatory affairs fax number: +41 32 566 54 09

Limitation of liability

Implantation of an intraocular lens is a surgical procedure that involves risks associated with eye surgery. SAV-IOL has provided the information and recommendations corresponding to this risk in addition to the requirements, methods and techniques used to implant the Lucidis Toric lens. SAV-IOL guarantees that the lens corresponds to the directions on the label for a period up to the expiry date. In particular, SAV-IOL accepts no liability in relation to any injury or damage suffered by the patient due to:

- The implantation method or technique used by the surgeon if the SAV-IOL recommendations are not followed.
- A lens prescription that is unsuitable for an individual patient.
- If the lens is used after the expiry date.

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