

EDEN

SWISS MADE INTRAOCULAR LENS

PREMIUM LENS FOR CATARACT SURGERY



EN



SWISS ADVANCED VISION
INTRAOCULAR LENS

A Unique Lens

Eden is best suited for patients with higher needs in near and intermediate vision such as reading, socializing, cooking, and watching television.

Enhancing the near and intermediate vision



About SAV-IOL

SAV-IOL is a Swiss company that designs and manufactures high-quality innovative intraocular lenses (IOLs) for cataract surgery.



SAV-IOL SA
Route des Falaises 74
2000 Neuchâtel
Switzerland



info@sav-iol.com



+41 32 566 54 00



sav-iol.com



SWISS ADVANCED VISION
INTRAOCULAR LENS

Important Information: SAV-IOL's lenses benefit from a lifetime of 20 years. The Instructions for Use (IFU) is available within the Eden box to provide detailed information about the lens. Cataract surgery is not completely risk-free. People with existing medical conditions are at a higher risk of developing complications. Other eye conditions before surgery may also impact the visual outcome. Inform your eye specialist immediately if you experience any unusual symptoms after the surgery.


What is a Cataract?

A cataract is a clouding of the eye's natural lens leading to a decrease in vision and potential blindness.

What's Next?


1

Assess the situation with a surgeon and select an artificial lens if surgery is needed.




2

Removal of the cataract through a small incision.



3

Placement of the artificial lens using an injector.




What is Eden?

Eden is a high-quality intraocular lens designed and made in Switzerland that replaces the eye's natural lens affected by a cataract.

Eden incorporates a unique EDOF optical technology offering patients an excellent continuous near to intermediate vision while minimizing loss of light and visual disturbances.

What is EDOF?

EDOF stands for Extended Depth of Focus. This is an optical design incorporated in intraocular lenses that aims to extend the depth of vision.



Benefits of Eden?

- ✓ Improves Light Intensity and Depth of Focus
- ✓ Reduces or Eliminates the Need for Glasses
- ✓ Enables a High-Quality Continuous Vision
- ✓ Optimizes the Near to Intermediate Vision

